

## How to upload BLOOD SUGAR METERS (EX. Accucheck, LifeScan/Verio, Ascensia/Contour or Freestyle)

1) Please visit this website: <https://diasend.com//ca-en> to create a Diasend account on your computer:

2) Once you have logged in, go to the tab "Account Admin" in your diasend® and add LMC's "Clinic ID" which is **33-15633**

Ask your health care provider for their Clinic ID and enter it here

The image shows a screenshot of the diasend account administration page. The page title is "Review my data" and the user name is "Jane Doe". The page has a navigation bar with "TOOLS", "ACCOUNT ADMIN", and "LOG OUT". Below the navigation bar, there are buttons for "Share data", "Connect app", "Recent Updates (Apr 2017)", and "Customer Support". The "Share data" section is active, showing options for "Share data", "Export data", "Subscription", and "Browse & connect apps". A message states: "By entering the Clinic ID of your health care provider team, they will have access to your medical data. Please consult your health care professional to receive their Clinic ID." Below this message, there is a "Data shared with: (24-83332)" and a "Stop sharing" button. The "Clinic ID" field is currently empty, with an "Update" button below it. To the right of the screenshot, there is a separate box labeled "Clinic ID" with an input field and a red "X" button. An orange arrow points from this box to the "Clinic ID" input field in the screenshot. Another orange arrow points from the "ACCOUNT ADMIN" tab in the screenshot to the "ACCOUNT ADMIN" label on the right.